

Zoom In and Out: A Comprehensive Immunologic Evaluation of Human, Murine and Rat Samples

Christos Nikolaou, Martin Lange, Alexandra Eichten, Oliver von Ahsen, Barbara Nicke, Charlotte Kopitz, Krzysztof Brzezinka

NUVISAN ICB GmbH | Muellerstr. 178 | 13353 Berlin | Germany


Introduction

Fluorochrome-conjugated antibodies applied in conventional FACS are widely used but have limited utility for high-parameter studies. Spectral cytometry overcomes those limitations since the emission spectrum of every fluorescence molecule is detected across a defined wavelength range. Employing spectral cytometry for immunophenotyping at NUVISAN we can zoom in and out of the immune system of:





1. Human
2. Mouse
3. Rat

with less sample material needed to extract complex information vital for immunological studies. Immune responses to cancer are highly influenced by the tumor microenvironment, where the delicate balance between suppressor and effector cells steers the prognosis after therapy. Proper establishment of multiparameter panels require accurate assessment of isolation and preservation protocols as practical considerations on epitope stability of immune cells vital for downstream assays and immunophenotyping.

We use spectral cytometry to map immune cell subsets and correlate them with tumor progression in human, mice and rats. Using unsupervised dimensionality reduction tools (e.g. viSNE) we identify major immune subsets as well as analyze their expression of stimulatory and inhibitory molecules in tissues and periphery.



Employing Spectral Cytometry for Immunophenotyping



Off the shelf NUVISAN panels available

Human Immunophenotyping Panel (≥30 parameters) of B, T, NK, DC, ILC and monocyte subsets in PBMCs.

Murine Immunophenotyping Panel (≥ 28 parameters) of B, T, NK, DC, and monocyte subsets in Thymus, blood, spleen, BM, tumors.

Rat Immunophenotyping Panel (≥ 13 parameters) of B, T, NK, DC, and monocyte subsets in Thymus, blood, spleen, BM, tumors.


human	mouse	rat
CD11c	CD4	CD4
CD45RA	CD19	CD45RA
CD3	Iy6G	CD8
CD25	CD69	CD3
IgD	CD45	CD62L
CD95	CD44	Gr
CD11b	CD11b	Ki67
CD38	CD62L	CD161a
CD57	Nkp46	CD45RC
CD27	CD137/41bb	Live/dead
CD123	CD8	CD28
CD127	CD45R/B220	IgM
HLADR	CTLA4	CD38
CCR7	F4/80	CD45
CD19	CD3	
CD16	Iy6C	
TCRgd	TCRgd	
CD14	CD11c	
CD8	PD1	
CD1c	CD25	
PD1		
CD56	SLAMF7	
CD45RA	CD206	
CD28	Live/dead	
SLAMF7	MHCII(A/E)	
CCR3	CD38	
CCR6		
cKIT		
IgM		
Live dead		

Advantage of Spectral flow cytometry:
Less sample material needed to extract more complex information!

Establishments & optimisations of FACS panels are crucial steps in flow cytometry experiments, as they optimize sensitivity, specificity, and reproducibility, leading to improved data quality and reliability. Find out more at <https://www.nuvisan.com/home.html> to see how we are committed to supporting cutting-edge research endeavours. As part of our comprehensive suite of services, we take great pride in offering specialized assistance as well as expertise in immunology and Fluorescence-Activated Cell Sorting (FACS) experiments.



Human immunophenotyping



Murine immunophenotyping

Experimental design 4 groups (aPD1 vs isotype control, aCTLA4 vs isotype control)

